

bioregen[®]

ULTRA PURE
HYALURONIC ACID

≈ REGENYAL
LABORATORIES

bioregen[®]

ULTRA PURE
HYALURONIC ACID

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The product is solely intended to be used by medical personnel. Do not use for applications other than indicated in this package leaflet.

DESCRIPTION

Bioregen is a resorbable medical device (sterile, apyrogenic and physiological gel) to be used for skin hydration and to correct superficial skin imperfections on face and body and in the dermal tissue repair process. The main component is Hyaluronic Acid sodium salt of non-animal origin, produced by bacterial fermentation.

COMPOSITION

Hyaluronic acid sodium salt 16 mg/g
Phosphate buffer, sodium chloride, water for injectable solutions... q.s. 1 g

PACKAGE 3x1 ml

- package leaflet
- 3 sealed blisters each containing 1 single-dose/single-use pre-filled sterile syringe
- adhesive tabs to be applied to the patient card for product traceability
- 3 needles

METHOD OF ACTION

Bioregen is injected into the dermis for the supplementation of the extracellular matrix and thanks to the hygroscopic properties of hyaluronic acid it restores the optimal physiological hydration of the skin.

INTENDED USE

Bioregen is a medical device produced in compliance with Directive 93/42/EEC and subsequent amendments and additions, intended for hydration of the skin and correction of superficial skin blemishes and in the dermal tissue repair process. It is particularly indicated in cases of scarring after superficial skin trauma (ex. acne and chicken pox scars) or stretch marks.

The outcomes of the procedure depend on skin type and the nature of the blemishes, with results being the better the less noticeable the nature of the defect.

DIRECTIONS FOR USE

Prior to any procedure with Bioregen, the doctor must conduct adequate anamnesis and an overall assessment of the patient's condition, to ensure the absolute absence of contraindications to the implant.

The areas to be treated must be identified and assessed based on criteria to do with lines and symmetry one must adhere to.

The procedure may be performed under local anaesthesia in order to ensure the necessary comfort to the patient.

The doctor must first inform the patient about the procedure, its nature, warnings, precautions, possible individual outcomes and potential adverse responses.

The area of the procedure should be cleaned with antiseptic solutions before implantation.

Extract the syringe from the blister, remove the top as shown in the figure and screw the needle or the cannula to the luer-lock, complete with protection.

Only remove the cover just before performing the procedure.

Bioregen is administered with a sterile needle compliant with Luer-Lock standards with standardised fittings, included in the package.

The injection should be performed in the top or middle dermis; the procedure is, however, at the physician's discretion and depends on the correction to be performed and the method implemented.

At the end of the procedure it is advisable to perform a gentle massage of the treated area for the product to be optimally distributed.

DOSE AND ADMINISTRATION

The volumetric markings on the syringes have indicative purposes, and the dosage to be used for the individual case is at the doctor's discretion.

The frequency with which to repeat the procedure depends on several factors, concerning both the patient's physiology (skin type, individual metabolism, anatomy, age) and lifestyle. Another factor to be taken into account is related to the chosen injection techniques. The procedure should be repeated periodically in order to maintain the results achieved, every 3-4 months.

PRECAUTIONS

The usual precautions in case of percutaneous procedures must be taken during the procedure.

The risks are those of common infection related to the type of treatment.

Bioregen must not be used on patients suffering from:

- infectious or inflammatory processes close to the area of the procedure
- known hypersensitivity to Keloids
- allergy to components
- immune system disorders
- chronic pathological skin condition
- conditions affecting the clotting factors or in the case of anticoagulant therapy in progress.

Around the time of treatment, the patient should avoid taking substances (Aspirin, NSAIDs, Vit. E) that act on blood fluidity, in order to reduce to the minimum the likelihood of bruising or bleeding of the injected areas.

The use of Bioregen is strictly excluded in areas where there are mammary, tendon, bone and muscle implants.

After treatment, and until swelling and redness go away completely, the treated areas should not be exposed to excessive heat (sun, UV tanning sessions, laser), or to intense cold.

After use dispose of needles and syringes according to the indicated procedures for hospital waste.

ADVERSE REACTIONS, WARNINGS AND CONTRAINDICATIONS

As is the case for percutaneous injections, there may be phenomena of:

- inflammatory reaction (redness, oedema, etc.), sometimes associated with itching and pain to the touch
- haematomas
- hardening or lumps at the injection site
- skin colouring or discolouring at the injection site

These phenomena generally resolve in few days. Should they persist more than a week, the patient must promptly contact his/her doctor.

Bioregen must be used intracutaneously and must not be injected into blood vessels.

No overdosage phenomena or interaction with medicinal products are known.

Do not use during pregnancy.

Bioregen is available in sterile single use packages.

It is forbidden to reuse the content for subsequent applications or on other patients.

The product must be used immediately after opening.

Any unused product must be disposed of.

Do not use the product if the package is damaged.

Do not mix with other injectables, nor use other implants in conjunction with Bioregen.

Perform the injection in a suitable facility observing the appropriate techniques.

Fill in the adhesive label included in the package and apply it to the patient card stored at the doctor's surgery and deliver the other to the patient. Keep away from the reach of children.

In rare cases, the product may cause an allergic reaction.

STORAGE METHODS

Bioregen must be stored between 2°C and 25°C.

Do not freeze. Do not expose to heat sources.

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REGENYAL LABORATORIES Srl